

**Leicester, Leicestershire and Rutland
Area Prescribing Committee**

LLR Continence Prescribing Formulary

October 2023

Date written/reviewed:	October 2023
Date of next review:	March 2025
Approved by LLR APC:	5.9.17
Version:	1.8
<i>If significant product price changes occur, product choice may be reviewed at an earlier date than indicated above.</i>	

Table of Contents

1.0 Introduction	3
2.0 Guidance for Prescribers.....	4
3. 0 Catheterisation Pack	6
4.0 Lubricating Gel for Catheterisation.....	7
5.0 Foley Catheters	8
5.1 Long term (up to twelve weeks)	8
5.2 Short / Medium Term Catheters – Up to 4 weeks use	10
6.0 Intermittent Catheterisation.....	11
6.1 Hydrophilic Intermittent Self Catheters.....	11
6.2 Pre-lubricated Hydrophilic Gel Intermittent Self Catheters	12
6.3 Discreet Pre-lubricated Hydrophilic Intermittent Self Catheters	13
6.4 Specialist Hydrophilic Intermittent Self Catheters with attached drainage bag	14
7.0 Catheter Valves	15
8.0 Catheter Drainage Bags	16
8.1 Leg Bags.....	16
8.2 Night Bags	17
9.0 Catheter & Leg Bag Accessories.....	18
9.1 Leg Bag Straps	18
9.2 Leg bag holders	18
9.3 Catheter Fixation Straps.....	18
9.4 Catheter Fixing Clips.....	19
10.0 Penile Sheaths	20
11.0 Catheter Maintenance Solutions	22
12.0 How to reduce waste	23
13.0 References	24
14.0 Acknowledgments.....	25

Printed copies may not be the up to date. Please use the electronic copies on
<https://www.areaprescribingcommitteeleicesterleicestershirerutland.nhs.uk/>

1.0 Introduction

The aim of this formulary is to identify the most appropriate use of urinary and intermittent catheters, penile sheaths, urine drainage bags and associated accessories for use in the community. It is expected that the majority of patients can be managed with the products recommended. This formulary provides guidance to prescribers for appropriate product choices for the majority of patients and is not intended to restrict patient choice.

Manufacturers' literature should be consulted for full details of indications of use, contra indications and side-effects.

A comprehensive continence assessment is required before considering any continence appliance; the emphasis should be on appropriate treatment.

Product selection should be made to meet patient needs on an individual basis as not all products are suitable.

Prices correct as of July 2017 Drug Tariff with the exception of penile sheaths which are correct of July 21 Drug Tariff.

Each section of the Continence Formulary will be reviewed by the designated Task and Finish Continence Team after three years unless significant circumstances arise that require immediate action. See individual sections for review dates.

2.0 Guidance for Prescribers

The table below provides a guide to prescribers on frequencies and quantities for continence appliances contained within this formulary. If you feel the patient is using any of the appliances excessively they can be referred to the Community Continence team for re-assessment.

Quantity guidance for continence appliances			
Appliance	Recommended frequency of change	Recommended <u>MONTHLY</u> quantity	Additional Information
Intermittent Catheters (IC)	Up to 5 times daily	Based on average use up to 150 per month (5 boxes of 30)	The frequency of use depends on the medical reason for the catheterisation, ranging from 1 to 5 times a day
Indwelling (Foley) Catheters	Long term (Up to 12 weeks)	One (1) every 3 months	The first time a patient is catheterised; two catheters should be prescribed (one for use and one for spare). Subsequent catheter prescriptions will be only for one catheter
	Medium Term (Up to 4 weeks)	One (1) every month	
Catheter Valves	Every 7 days	Five (5)	Pack of 5 per month
Leg Bags	Usually one bag every 5 – 7 days	Five (5) bags	1 box of 10 every 2 months
Night Bags (Single Use)	One bag per day	Thirty (30) bags	-
Night Bags (Drainable)	Usually one bag every 5 – 7 days	Five (5) bags	1 box of 10 every 2 months
Leg bag straps	-	One (1) pair supplied with leg bags as standard	Pack of 5 pairs every 6 – 9 months
Leg bag holders	-	-	Pack of 4 every 6 months
Catheter fixation clips	Usually every 7 days	Five (5)	Pack of 5 per month
Catheter retaining straps	-	-	Pack of 5 every 4-6 months
Penile Sheaths	Usually one per day	Thirty (30)	-

Patients requesting continence appliances and accessories

Prescribers should be aware that patients receive promotional material and samples of many continence appliances and accessories via magazines and conventions. Requests for such items should only be considered if recommended by the continence nurse.

Dispensing of prescriptions

It is important that patients are aware they have a choice as to which dispensing appliance contractor (DAC) or community pharmacy dispenses their appliance prescription and that they can change their choice at any time.

Retrospective prescriptions

Retrospective prescriptions should not be issued by a prescriber.

In the case of an emergency, the DAC or pharmacy contractor may supply an appliance only at the request of a prescriber. The prescriber should undertake to give the DAC or pharmacist a non-electronic prescription form or non-electronic repeatable prescription in respect of the appliance requested within 72 hours. Alternatively, the prescriber must transmit to the ETP service an electronic prescription within 72 hours.

If a contractor requests a prescription “post supply” GPs are entitled to refuse the supply of a prescription

3.0 Catheterisation Pack

Last reviewed: April 2024

Next review: April 2027

Cath-It Catheter Insertion/Removal Pack S/M/L (catheter, leg bag, sodium chloride and lubricating gel not included)- For use in CLINIC SETTING

Product	Manufacturer Code	Pack Size
Cath-It Catheter Insertion/Removal Pack (Small Gloves)	908410 (EHC522)	Pack of 5
Cath-It Catheter Insertion/Removal Pack (Medium Gloves)	908420 (EHC523)	Pack of 5
Cath-It Catheter Insertion/Removal Pack (Large Gloves)	908430 (EHC524)	Pack of 5

Prosys® Vesica Community catheterisation Pack (catheter and lubricating gel not included)- For use in COMMUNITY SETTING

Product	Manufacturer Code	Pack Size
Catheterisation pack (slide tap leg bag)	Vescp-01	1
Catheterisation pack (lever tap leg bag)	Vescp-02	1
Catheterisation pack (slide tap leg bag)	Vescp-03	1

- Facilitates aseptic non touch technique (ANTT) to reduce risk of catheter associated urinary tract infections (CAUTI). Additional material present in a clinic setting will facilitate this when utilising the ***Cath-It Catheter Insertion/Removal pack***.
- Promotes and standardises best practice using a two layer system, layer 1 catheter removal kit, layer 2 catheter insertion kit.
- Cost-effective catheterisation pack as supplied in one packet reducing nursing time and FP10 costs.

4.0 Lubricating Gel for Catheterisation

Next review: October 2023

- These should be prescribed by a registered clinician and authorisation recorded in patients notes.
- Provides a protective layer to ensure smooth passage of a catheter, reducing potential injury to the urothelium and subsequent risk of possible urethral damage.
- It can facilitate pain free insertion.
- Helps to reduce the risk of associated infection.
- One to be used at each catheter change.
- For urethral catheter insertion: 6mls female, 11mls male.
- Licensed for suprapubic use where appropriate: 6ml to be used.
- Be aware of the potential for an anaphylactic reaction to chlorhexidine

DO NOT USE in patients with a known allergy to chlorhexidine (MHRA Warning, 2012)

Manufacturer	Product	Order Code	Pack Size
CliniMed	Instillagel® (Lidocaine 2% & Chlorhexidine 0.25% gel)	6ml (Female)	1
		11ml (Male)	1

5.0 Foley Catheters

Next review: October 2023

- Patient assessment, monitoring and ongoing support is essential to prevent urinary tract infections and improve quality of life for those patients who have a urinary catheter.
- All patients should have a **catheter passport** which is completed and kept up to date.
- Aim to remove catheter as soon as possible following insertion, and review continuing need. Also aim to use the smallest size that provides adequate drainage to avoid problems such as bypassing.
- Select the correct length and type of catheter; Standard / females & also suitable for suprapubic catheterisation. Female - **for females only**
- Use the correct amount of sterile water to fully inflate the balloon as per manufacturer's instructions (usually 10ml). Under-inflation can distort the angle of the catheter tip, causing bladder spasm. The trusts policy states only catheters with 10ml balloons should be used.

Choice of product depends upon:

1. Assessment & diagnosis
2. Patient choice
3. Local guidelines

If problems persist i.e. frequent re-catheterisations, bypassing, blocking, **contact the community continence service for advice/support.**

Intermittent self-catheterisation is considered the gold standard treatment and should be used in preference to an indwelling catheter if it is clinically appropriate and is a practical option for the patient. Please refer these patients to the continence service for individual assessment where appropriate.

A prescription must be generated and authorised by the prescriber prior to any ordering/delivery of catheters, this will prevent excessive quantity and frequency of ordering

5.1 Long term (up to twelve weeks)

- Packaged in single units.
- Prescribe two initially and from then on only one to be prescribed at a time.
- Catheters listed are licensed for both urethral/suprapubic uses.

If catheter life is less than four weeks i.e. requiring recurrent re-catheterisation, consider a medium term catheter.

Teleflex Rüsch Brilliant AquaFlate®: All-Silicone with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation

Product	Order Code	Balloon Size (ml)	Gauge (Ch)	Pack Size
STANDARD / MALE	DA310112	10	12	1
	DA310114		14	
	DA310116		16	
	DA310118		18	
	DA310120		20	
	DA310122		22	
	DA310124		24	
FEMALE	DA210112	10	12	1
	DA210114		14	
	DA210116		16	
	DA210118		18	
	DA210120		20	
	DA210122		22	
	DA210124		24	

Coloplast Folysil® X-Tra: Silicone open ended catheter with pre-filled syringe for balloon inflation and empty syringe for balloon deflation. **This catheter should be reserved and trialled for patients whose catheter blocks frequently.**

Product	Order Code	Balloon Size (ml)	Gauge (Ch)	Pack Size
MALE	AA8C12	10	12	1
	AA8C14		14	
	AA8C16		16	
	AA8C18		18	

5.2 Short / Medium Term Catheters – Up to 4 weeks use

Teleflx Rusch PTFE Aquaflate®: PTFE Coated latex with sterile water filled syringe for balloon inflation and empty syringe for balloon deflation

Product	Order Code	Balloon Size (ml)	Gauge (Ch)	Pack Size
STANDARD / MALE	DP310112	10	12	1
	DP310114		14	
	DP310116		16	
	DP310118		18	
	DP310120		20	
	DP310122		22	
	DP310124		24	
FEMALE	DP2101	10	12	1
	DP2101		14	
	DP2101		16	
	DP2101		18	
	DP2101		20	
	DP2101		22	
	DP2101		24	

6.0 Intermittent Catheterisation

Next review: October 2023

Intermittent catheterisation (IC) is a technique that may be carried out by the individual, their carer (s) or by appropriately trained healthcare staff. IC has a lower risk of infection compared with indwelling catheterisation, and so should be considered in preference to an indwelling catheter (NICE 2010)

Prior to implementing intermittent catheterisation you need to ensure that an assessment has been completed and that:

- The bladder has the capacity to store urine adequately between catheterisations
- Individuals and/or carer (s) are motivated and have the manual dexterity and the ability to carry out the procedure.

- The frequency of change for intermittent catheters can vary by patient. On average these should be used up to five times a day.
- Patients who are catheterising more often by this method may need a review of suitability and training technique reinforced.

6.1 Hydrophilic Intermittent Self Catheters (FIRST LINE)

These should be used first line in patients self-catheterising and either contains a sterile water sachet within the packaging, which is required to activate the catheter prior to use or is pre-lubricated in a hydrophilic coating.

The patient should be shown the catheters and given the choice of which they would prefer to use.

Male: Hydrophilic ISC in saline solution

Product	Manufacturer Code	Gauge / Ch	Pack Size
Coloplast Speedicath®	28408	8	30
	28410	10	
	28412	12	
	28414	14	
	28416	16	
	28418	18	
WellSpect Healthcare LoFric® Origo™ (40cm)	4431025	10	30
	4431225	12	
	4431425	14	
	4431625	16	
	4431825	18	

Female: Hydrophilic ISC in saline solution

Product	Manufacturer Code	Gauge / Ch	Pack Size
Coloplast Speedicath®	28506	6	30
	28508	8	
	28510	10	
	28512	12	
	28514	14	
	28516	16	
WellSpect Healthcare LoFric® Sense™ (15cm)	4160825	8	30
	4161025	10	
	4161225	12	
	4161425	14	

6.2 Pre-lubricated Hydrophilic Gel Intermittent Self Catheters (RESTRICTED)

These catheters should be reserved for patients who are unable to break the water sachet and require a pre-lubricated catheter to self-catheterise. These catheters are coated in glycerine / water mix to aid insertion.

Male: Pre-lubricated hydrophilic gel ISC: Actreen® Hi-Lite (B Braun Medical)

Product	Manufacturer Code	Gauge / Ch	Pack Size
Actreen® Hi-Lite (41cm)	238208E	8	30
	238210E	10	
	238212E	12	
	238214E	14	
	238216E	16	

Female: Pre-lubricated hydrophilic gel ISC: Actreen® Hi-Lite (B Braun Medical)

Product	Manufacturer Code	Gauge / Ch	Pack Size
Actreen® Hi-Lite (20cm)	238306E	6	30
	238308E	8	
	238310E	10	
	238312E	12	

	238314E	14	
	238316E	16	

6.3 Discreet Pre-lubricated Hydrophilic Intermittent Self Catheters (RESTRICTED)

These catheters are suitable for patients who require a discreet pre-lubricated catheter as they are unable to break open the water sachet contained in first line choices. Please note both the LoFric® Origo™ and LoFric® Sense™ catheters are discreet catheters that can be folded, however these require the water sachet to be broken prior to use.

Discreet Male: Pre-lubricated hydrophilic ISC: Speedicath® Flex (Coloplast Ltd)

Product	Manufacturer Code	Gauge / Ch	Pack Size
Speedicath® Flex	28920	10	30
	28922	12	
	28924	14	
	28916	16	

Discreet Female: Pre-lubricated hydrophilic ISC: Speedicath® Eve (Coloplast Ltd)

Product	Manufacturer Code	Gauge / Ch	Pack Size
Speedicath® Compact Eve	28110	10	30
	28112	12	
	28114	14	

6.4 Specialist Hydrophilic Intermittent Self Catheters with attached drainage bag **(RESTRICTED)**

These catheters are hydrophilic and contain a water sachet for activation prior to use, they also have an attached urine drainage bag. These catheters should be reserved for patients who are unable to transfer to a toilet and/or being catheterised by another individual. If the patient / carer are unable to break the water sachet the Continence Team will assess the product that is best suited for this individual.

Product	Manufacturer Code	Gauge / Ch	Pack Size
LoFric® Hydro-Kit™ <i>(MALE 40cm)</i>	4200825	8	20
	4201025	10	
	4201225	12	
	4201425	14	
	4201625	16	
	4201825	18	
LoFric® Hydro-Kit™ <i>(FEMALE 20cm)</i>	4230825	8	20
	4231025	10	
	4231225	12	
	4231425	14	
	4231625	16	
	4231825	18	

7.0 Catheter Valves

Last reviewed: Oct 2023

Next review: Oct 2026

- For use with indwelling catheters only.
- These provide a discrete alternative to drainage bags.
- Their use helps to imitate normal bladder function by allowing the bladder to fill and empty, maintaining normal capacity and tone.
- They allow the catheter balloon to be lifted from the bladder wall decreasing the risk of bladder wall erosion and trauma to the bladder neck. (*RCN guidelines, 2012*)
- To be considered 2-3 weeks prior to trial without catheter to regain bladder function and tone.
- Catheter valves should be changed every 5-7 days.
- No more than one packet (5) should be prescribed every month.

Contraindications:

- Reduced bladder capacity
- No bladder sensation
- Cognitive impairment
- Poor manual dexterity
- Renal impairment
- Post Radical Prostatectomy

Manufacturer	Product	Order Code	Pack Size
Clinisupplies Ltd	Prosys® Catheter Valve	PCV3942	5

8.0 Catheter Drainage Bags

Next review: October 2023

8.1 Leg Bags

- Leg bags may be worn in different positions on the leg i.e. inside thigh/calf, which is an individual choice and which will determine the length of the inlet tube.
- Leg bags should be changed every 5-7 days. (*Manufacturer's recommendation as detailed in DoH Drug Tariff July 2017*)
- Use ANTT when changing a leg bag advocated by Trust.
- Maintaining a closed drainage system (i.e. not removing the leg bag when attaching a night bag) reduces the risk of infection.
- The leg bag must remain connected to the catheter and linked to the night bag if additional drainage capacity is required overnight.
- It is advisable to alternate leg each day to prevent sores.
- 500ml bags have been listed below, as this is the most common size used. **Other volume leg bags are available.**

Prosyst[®] Sterile Leg Bags with lever tap: Containing 1 pair of non-latex gloves and soft elasticated cotton velcro straps

Manufacturer	Product	Order Code	Pack Size
CliniSupplies LTD	Prosyst [®] Sterile Leg bag Short tube with Lever Tap	P500S-LT	10
	Prosyst [®] Sterile Leg bag Long tube with Lever Tap	P500L-LT	

Prosyst[®] Sterile Leg Bags with slide action tap: Containing 1 pair of non-latex gloves and soft elasticated cotton Velcro straps

Manufacturer	Product	Order Code	Pack Size
CliniSupplies LTD	Prosyst [®] Sterile Leg bag Short tube with Slide action tap	P500S	10
	Prosyst [®] Sterile Leg bag Long tube with Slide action tap	P500L	

8.2 Night Bags

Non Drainable Night Bags (single use, non-sterile)

- Suitable for night time use for the collection of urine from indwelling catheters.
- These should be used in all community and in patients.
- The position of the bag should be below bladder level to enhance drainage.
- Night bags should be directly connected to the leg bag to maintain a closed system. They should not be used directly connected to the catheter.
- Night bags should be used in conjunction with a catheter stand, available FOC from CliniSupplies. (not available on FP10) Please call 0800 0854957 to obtain one.
- Nursing/residential homes should always use a single use, non-drainable night bag attached to a leg bag.
- Non drainable night bags are preferred option as single use reduces the risk of infection.
- No more than 3 boxes of 10 should be issued per month.

Manufacturer	Product	Order Code	Pack Size
CliniSupplies Ltd	Prosys® Non-sterile 2 litre night bag with single use drainable tap	PSU2	10

Drainable Night Bags (Re-usable)

- For community bed bound patients it may be appropriate for a sterile drainable 2 litre bag to be connected directly to the catheter. (Not recommended for nursing/residential home patient due to cross infection)
- Drainable night bags should be changed every 5-7 days. (*Manufacturer's recommendation, DoH Drug tariff July 2017*)
- A box of 10 drainable night bags should last two months.

Manufacturer	Product	Order Code	Pack Size
CliniSupplies Ltd	Prosys® Sterile 2 litre night bag with slide outlet tap	P2000	10
CliniSupplies Ltd	Prosys® Sterile 2 litre night bag with lever outlet tap	P2000-LT	10

9.0 Catheter & Leg Bag Accessories

Next review: October 2023

9.1 Leg Bag Straps

- It is extremely important that both the catheter and leg bag is well supported to reduce traction and trauma to the bladder neck/urethra.
- There is 1 pair of CliniSupplies Prosys® leg bag straps provided in each box of leg bags. Should the patient require more these can be ordered separately.
- A pack of 5 pairs of leg bag straps should last six to nine months.

Manufacturer	Product	Order Code	Pack Size
CliniSupplies LTD	Prosys® soft elasticated cotton leg bag support straps with silicone leg grip™	P10LS	10 (5 pairs)

9.2 Leg bag holders

- Ugo Fix Sleeves can be used as an alternative to, or alongside, leg straps. They are particularly good if patients have frail skin or problems with straps digging into or rubbing against legs, as they distribute the weight of the bags more uniformly.
- A pack of 4 Ugo Fix sleeves should last four to six months.

Manufacturer	Product	Order Code	Pack Size
Optimum Medical Solutions	Ugo Fix sleeve leg bag holder	Small (24-39cm) - 3005	4
		Medium (36-55cm) – 3006	
		Large (40-70cm) - 3007	
		X-Large (65-90cm) - 3008	

9.3 Catheter Fixation Straps

- The Simpla® G-strap is used as a retaining strap which secures the tubing or catheter firmly and comfortably against the leg acting as a shock absorber for all indwelling catheters.
- A pack of 5 Simpla® G straps should last four to six months.

Manufacturer	Product	Order Code	Pack Size
--------------	---------	------------	-----------

Coloplast Ltd	Simpla® G Strap	383001	5
---------------	-----------------	--------	---

9.4 Catheter Fixing Clips

- The Ugo fix gentle is an alternative to the Simpla® g-strap to ensure the catheter or tubing is securely safely.
- These are designed to be moved and placed multiple times. The adhesive allows one clip to remain in place for up to 7 days.
- A pack of 5 Ugo fix gentle catheter clips should last at least one month.

Manufacturer	Product	Order Code	Pack Size
Optimum Medical Solutions	Ugo Fix Gentle	3004	5

10.0 Penile Sheaths

Last reviewed: Sept 2021

Next review: Sept 2024

- Can offer a valuable alternative method of urinary incontinence management for men.
- It is important to accurately assess for type and size of sheath using manufacturers measuring device.
- Sheaths can be worn for up to 24hrs.
- Over ordering more than one box of thirty per month may indicate poor fit - refer to Community Continence Service for advice.

First line penile sheaths

Place in Formulary	Manufacturer	Product	Order Code	Size (mm)	Pack Size
First line	Clinisupplies Ltd	Prosys flofit Standard	24	PSS24STD	30
			28	PSS28STD	
			31	PSS31STD	
			35	PSS35STD	
			40	PSS40STD	
First Line	Clinisupplies Ltd	Prosys flofit Pop on	24	PSS24POPON	30
			28	PSS28POPON	
			31	PSS31POPON	
			35	PSS35POPON	
			40	PSS40POPON	
First line	Clinisupplies Ltd	Prosys flofit Wideband	24	PSS24WIDEBD	30
			28	PSS28WIDEBD	
			31	PSS31WIDEBD	
			35	PSS35WIDEBD	
			40	PSS40WIDEBD	

Second line penile sheaths

Place in Formulary	Manufacturer	Product	Order Code	Size (mm)	Pack Size
Second Line	Greatbear Healthcare Ltd	GB libra Standard	24	GBLSS24	30
			28	GBLSS28	
			31	GBLSS31	
			33	GBLSS33	
			35	GBLSS35	
			40	GBLSS40	

Second line	Greatbear Healthcare Ltd	GB libra Pop on	24	GBLSP24	30
			28	GBLSP28	
			31	GBLSP31	
			33	GBLSP33	
			35	GBLSP35	
			40	GBLSP40	
Second Line	Greatbear Healthcare Ltd	GB libra Wideband	24	GBLSW24	30
			28	GBLSW28	
			31	GBLSW31	
			33	GBLSW33	
			35	GBLSW35	
			40	GBLSW40	

11.0 Catheter Maintenance Solutions

Next review: October 2023

- Catheter maintenance solutions should be used as part of a treatment plan and **not** as routine for every patient.
- The solution should only be used following a patient assessment by a Registered Health Care Professional who will then arrange an authorisation and prescription.
- Solutions should only be considered for short term use to prevent encrustation in indwelling catheters or to dissolve crystal formation prior to the removal of a catheter to prevent urethral trauma.
- They are treatment preparations for dissolving encrustation only and should only be prescribed for 4 weeks, then reviewed.
- **They are not to be used routinely or to attempt to unblock a non-draining catheter.**
- It is good practice to cut open the catheters on removal, from those patients where blocking is problematic to see if the lumen is blocked by sediment deposit.
- The total amount of material dissolved when two solutions are used sequentially is significantly better than using just a single solution. *(Getliffe 2000)*

Types of solutions

- **Sodium chloride (S)** Sodium chloride may be helpful for removing blood clots, mucus and any other mechanical blockages.
- **Citric acid 3.23% (G)** will dissolve encrustation. It is recommended that two instillations should be used sequentially. The agitation property of this product should be used with care.
- **Citric acid 6% (R)** will dissolve crystals and is twice the strength of the G solutions. It is recommended that two instillations should be used sequentially. The agitation property of this product should be used with care.

Product	Order Code	Volume	Pack Size
Optiflo® S (0.9% saline)	CSS50	50ml	1
Optiflo® G (3.23% Citric Acid Suby G)	CSG50	50ml	1
Optiflo R® (6.0% Citric Acid Solution R)	CSR50	50ml	1

12.0 How to reduce waste

Next review: October 2023

- Any continence accessories (i.e. leg bag straps, leg bag holders and catheter retaining straps) should be prescribed as an **ACUTE** to avoid continual re-ordering.
- Accessories can be washed and re-used. Refer to manufacturer's guidance.
- Long term (12 week) foley catheters should be prescribed on a 3 monthly basis after initial prescription.
- Legs bags should only be requested and issued on a 2 monthly basis as 1 box of 10 lasts 2 months.
- Use quantity guidance (pg. 4) as a guide if patients are ordering more frequently.
- Do not supply prescriptions of new items unless requested or authorised by the continence nursing team.
- Ensure that only one contractor is requesting appliances. They can be supplied from either a DAC or a Pharmacy **NOT BOTH**.
- Prior to request the contractor should check with the patient whether an item is required. If you feel a patient may be over ordering do a spot check on their stock.
- Patients should receive dry wipes and disposable bags with each delivery FOC – **these should not be prescribed**.
- If a patient is ordering more than the recommended quantity of intermittent catheters or penile sheath this may indicate poor technique or fit. Speak to the Community Continence team for support.

13.0 References

1. Department of Health - Drug Tariff July 2017
2. PrescQIPP (<https://www.prescipp.info/>)
3. Getliffe K, Hughes S and Claire M (2000) The dissolution of urinary catheter encrustation British Journal Of Urology 85(1) .: 60-4
4. Guide To Urine Drainage Systems - Coloplast - 2010
5. NICE Clinical Guidelines 139 (2012) Prevention and Control of Healthcare associated infections in primary and community care
6. NICE Clinical Guidelines 171 (2013) - Urinary incontinence in women - management
7. RCN Guidelines (2012) Catheter care –RCN guidance for Nurses
8. Turner B and Dickens N (2011) Long Term Urethral catheterisation Nursing Times Vol 21 no 4 32-38
9. Worcestershire Continence Care Formulary (2017)
10. MHRA Drug alert (Oct 2012) All medical devices and medicinal products containing chlorhexidine - Risk of anaphylactic reaction due to chlorhexidine allergy for information
11. Staffordshire Continence Prescribing Formulary 2017
12. National Patient Safety Agency (2009) Rapid response report: Female urinary catheters causing trauma to adult males
13. Roberts JA, Fussell EN, Kaack MB (1990) Bacterial adherence to urethral catheters. Journal of Urology. 144, 2 Pt 1, 264-269.
14. Loveday HP, Wilson JA, Pratt RJ, Golsorkhi M, Tingle A, Bak A, Browne J, Prieto J, M. Wilcox (2014) epic3: National Evidence-Based Guidelines for Preventing Healthcare-Associated Infections in NHS Hospitals in England. Journal of Hospital Infection. 86S1, S1 – S70.

14.0 Acknowledgments

Special thanks for participation in development of the formulary to;

Vishal Mashru	Senior Pharmacist Medicines Optimisation West Leicestershire CCG
John Fletcher	Independent Nurse Advisor
Rosalind Gretton	Continence Lead Nurse (LPT)
Chris Rippin	Continence Lead Specialist Continence Service (LPT)
Janet Browning	Specialist Nurse (UHL)
Kash Bhayani	Patient Representative
Claire Bennett	Deputy Sister Urology (UHL)
Mini Anto	Urology Nurse Specialist (UHL)
Jaskarn Rai	Urology Consultant (UHL)
Elinor Howcroft	Continence Lead Nurse Specialist (UHL) Penile Sheath update
Cathy Oakley	Continence Nurse (LPT) Penile Sheath update
Sarah Coombes	Continence Specialist Nurse (UHL) Penile Sheath update
Stuart Punnett	Formulary Pharmacist (LLR) Penile Sheath update
Rebecca Dickenson	Healthcare Community Interface Pharmacist (LLR) Penile Sheath update
Vinod Patel	Specialist Medicines Information Technician (UHL) Penile Sheath update

Version	Section	Description of amendments	Date	Author/amended by
1.0	All	New guidance	10/8/17	Vishal Mashru
1.1	All	Final version post comments from formulary group	25/8/17	Vishal Mashru
1.2	2.0	Medium term catheter recommended monthly quantity amended to 1 every month	12/4/18	Vishal Mashru
1.3	3.0 & 6.4	Codes updated	17/07/18	Vishal Mashru
1.4	10.0	Formulary choices for penile sheaths updated	13/9/21	Rebecca Dickenson
1.5	7.0	Interim catheter valve added due to discontinuation of the previous valve	10/10/22	Rebecca Dickenson
1.6	7.0	Catheter valve after formulary work	3/10/23	Rebecca Dickenson/Cathy Oakley
1.7	3.0	Formulary review of catheterisation packs	04/04/24	Nazia Suleman/Cathy Oackley/Chris

				Rippin/ Jen Zulfugarova
1.8	All	Extension to review date to March 2025	21/10/24	Amandeep Hayer